

DPM1B

Pneumatic Transducer Tester

Technical Data



The DPM1B pneumatic transducer tester is designed to measure the positive and negative pressures of medical devices in either liquid or gaseous form, and to generate pressure within the ± 300 mmHg range to assist in repair and quality control.

Key features

- Battery operated
- Generates and measures positive or negative pressures
- Operates with gas and liquid
- Troubleshooting with 1 % accuracy

Technical specifications

Pressure measurement

Operating range: ± 300 mmHg

Pressure generation range: ± 300 mmHg

Accuracy: ± 1 % of reading or ± 1 mmHg

Units of measure: mmHg

Environmental requirements

Operating temperature: 10 °C to + 40 °C

(+ 50 °F to + 104 °F)

General information

Display/Control: 0.5 in LCD with LO BATT indication

Power: 9 V alkaline battery

Weight: 10 oz (260 g)

Ordering information

Item numbers/Descriptions

2249779 DPM 1B Pneumatic Transducer Tester

Standard accessories

2572314 User/service manual

2242666 4-way stop cock adapter

614487 9 V alkaline battery

Optional accessories

2242666 4-way stop cock adapter

2212678 Acrodisc filter

About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance.

Highly credentialed and equipped with a NVLAP Lab Code 200566-6 accredited laboratory, Fluke Biomedical also offers the best in quality and customer service for all your equipment calibration needs.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical

6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe

Science Park Eindhoven 5110
5692EC Son, The Netherlands

For more information, contact us:

In the U.S.A. (800) 850-4608 or
Fax (440) 349-2307

In Europe/M-East/Africa +31 40 267 5435 or
Fax +31 40 267 5436

From other countries +1 (440) 248-9300 or
Fax +1 (440) 349-2307

Email: sales@flukebiomedical.com

Web access: www.flukebiomedical.com

©2006-2009 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A.
11/2009 2817969B D-EN-N

Modification of this document is not permitted without written permission from Fluke Corporation.